

Accessories for Pet Injections & Blood Sample Collection



Inject-Ease™

Pet Injections Made Easy

Re-usable device that makes pet subcutaneous injections easy. Simply place a loaded syringe in the Inject-Ease, place the tip against the skin, and press the button to automatically deliver the needle through the skin. You control the rate at which the medicine is injected.

- Usable with subcutaneous injection treatments, including diabetes and other conditions
- Special tip & technique- helps reduce the pain by masking the needle puncture
- One button use provides aid to the pet owner with limited dexterity and coordination
- One press of the button inserts the needle through the skin automatically and easily
- Allows use of different size insulin syringes-works with 30, 50, or 100 unit syringes
- Spacer rings allow to vary the injection depth
- Quiet so it does not startle the pet
- Easy to practice injection technique to increase confidence
- Stabilizes the device on the skin to reduce early withdrawal and loss of fluid
- Affordable with 2 year warranty
- Part# IE400V




Insul-eze™

Syringe Magnifier

One-piece syringe magnifier that holds syringe and insulin bottles in place for those with visual or manual dexterity problems.

- Aids in syringe loading for visually impaired
- Magnifies all syringe calibrations two times for easy reading
- Holds syringe and bottle firmly while drawing dosage
- One piece unit accommodates most syringes and bottles
- Can be attached to refrigerator or wall with Velcro or magnet for easy viewing and loading
- Part #IZ6000




Insul-Cap™

Stores, Identifies Insulin Bottles, and Helps Load Syringes

Insulin bottle caps designed to aid in storing, identifying and making drawing insulin into a syringe easier.

- Allows one-handed syringe loading
- Fits most standard syringes and insulin bottles
- Easy to put on, store and remove
- Keeps bottle tops sanitary
- Safely holds insulin bottle while user holds syringe
- Part #IC200S



All products are trademark of
AMBI, INC. All rights reserved.

AMBIVETPRODUCTS
World Class Veterinary Products

908 Sir Francis Drive
Capitola, CA 95010, USA
Phone/Fax: 831-475-1765
email: www.ambimedinc.com
web: info@ambimedinc.com

Specifications are subject to change with-
out notice. ©2010 AMBIVETPRODUCTS.
All rights reserved.

Insul-Tray™ *Insulin Bottle Tray*

Small and sturdy Insulin bottle holder that keeps bottles organized and standing upright.

- Fits most standard insulin bottles
- Easy to put on, store and remove
- Part #ITR500S

EZ-Lets™ *Lancets*

High-quality lancets compatible with most standard lancing devices.

- Precision-ground tri-bevel point for clean, precise punctures
- Full range of sizes available — 21 gauge, 26 gauge, and 28 gauge

auto-Lancet *Adjustable Tip Lancing Device*

Re-usable quiet and comfortable lancing device that assists in obtaining capillary blood samples for testing.

Five level adjustable tip allows variation of lancet penetration

- Pen shape available in standard (5") or mini (3.5") size
- Works with industry standard lancets
- Five year warranty

Ordering Information

AMBIVETPRODUCTS Order Hotline: 831-475-1765 E-mail: info@ambivetproducts.com

Please note that not all products are available in all markets. For full details please contact your local reseller or AMBIVETPRODUCTS.

Special Orders

For special configurations and high volume orders, contact AMBIVETPRODUCTS or send e-mail to info@ambivetproducts.com.