

EZ Lets™

Premium Blood Lancets



Experience the best in blood lancets from AMBIMEDINC, a global leader in the design and manufacture of medical device products focused on reducing pain. EZ-Lets are precision engineered for the finest in needle edges, reducing discomfort and encouraging the healing process. EZ-Lets fit most re-usable finger lancing devices, and are available in four sizes: 21 gauge, 26 gauge, 28 gauge and 30 gauge.

Rely on AMBIMEDINC to deliver the highest quality blood lancets for capillary sampling requirements!

Key Features

Safety

- Sterile blood lancets
- Single use to eliminate cross-contamination
- Twist cap to expose needle
- Cap may be used to receive used needle, to avoid blood contamination

Precision Engineered

- Finely-ground needle edges for cleaner lance
- Fine points for least amount of capillary damage
- Consistent quality

Key Benefits

- Four sizes available: 21 gauge, 26 gauge, 28 gauge and 30 gauge.
- Available in 100 count and 200 count boxes
- Ultimate ease of use — insert in device and remove needle cap to use
- Highest quality for maximum comfort and least amount of cell damage

Ordering Information

AMBIMEDINC Order Hotline:
1-831-475-1765
E-mail: info@ambimedinc.com

Please note that not all products are available in all markets. For full details, please contact your local reseller or AMBIMEDINC

Special Orders

For special configurations and high volume orders, contact AMBIMEDINC or send e-mail to info@ambimedinc.com

EZ Lets is a trademark of AMBIMEDINC.
All rights reserved.

*Other trademarks are property of their respective holders.

Made in the U.S.A.

AMBIMEDINC
World Class Healthcare Products

908 Sir Francis Avenue
Capitola, CA • 95010, USA
Phone/Fax 1-831-475-1765
e-mail: info@ambimedinc.com
web: www.ambimedinc.com

Specifications are subject to change without notice.
© 2015, AMBIMEDINC. All rights reserved.

Specifications & Packaging Configurations

30 Gauge

Needle size: 30 gauge
Point configuration: Tri-bevel

30 gauge

PN/REF: ETP100 (100 ct.)
UPC Code: 0 25294 00825 9

PN/REF: ETP200 (200 ct.)
UPC Code: TBD

28 Gauge

Needle size: 28 gauge
Point configuration: Tri-bevel

28 gauge

PN/REF: ETU100 (100 ct.)
UPC Code: 0 25294 00820 4

PN/REF: ETU200 (200 ct.)
UPC Code: TBD

26 Gauge

Needle size: 26 gauge
Point configuration: Tri-bevel

26 gauge

PN/REF: ETY100 (100 ct.)
UPC Code: 025294 00815 0

PN/REF: ETY200 (200 ct.)
UPC Code: TBD

21 Gauge

Needle size: 21 gauge
Point configuration: Tri-bevel

21 gauge

PN/REF: ETO100 (100 ct.)
UPC Code: 0 25294 00810 5

PN/REF: ETO200 (200 ct.)
UPC Code: TBD

All Models

Length (w/tab): 1.5 in. (38.1 mm)

Diameter: 0.25 in. (6.35 mm)

Materials: Polypropylene (needle carrier), stainless steel (needle), Latex free

For use with most industry-standard finger-lancing devices, including:

auto-Lancet™, SoftTouch™, Qwik-Let™, Penlet® Plus, Penlet® II,
MicroLet™ Vaculance, Dialet™, OneTouch UltraSoft™, Monojector™, Freestyle®,
B-D • Medisense®, Microlet™, Autolet®